

HALOGEN-LICHTQUELLE DI LITE

Hellstes Licht für viele Anwendungsbereiche
*Efficient Light Source for Use in a Lot of
Applications*





HALOGEN-LICHTQUELLE DI LITE Art.-Nr. 97.300.00

- Hochleistungslichtquelle mit Halogenlampen für maximale Ausleuchtung
- Einfache Bedienung
- Leistungsstark für höchste Anforderungen in der Medizin
- Mit 2 Lichtentnahmestellen
- Optimal geeignet in Verbindung mit 1 CCD Kameras
- Inklusive Komplettzubehör

Technische Daten:

Lampen	2 x Halogen 250 W
Farbtemperatur	ca. 3200 K
Lichtintensitätsregelung	mechanische Blende (stufenlos)
Lichtleiteranschluß	System Storz
Netzspannung	115/230/240 V AC
Frequenz	50/60 Hz
Leistungsaufnahme	260 W
Betriebstemperatur	+10°C bis +40°C
Maße (BxHxT)	300 x 145 x 320 mm
Gewicht	ca. 10 kg

Technische Änderungen vorbehalten.

- *High performance light source with 2 halogen lamps for optimized light output*
- *Easy to operate*
- *Meets the highest light requirements for operative procedures*
- *2 light guide ports*
- *Recommended for use with 1 Chip CCD cameras*
- *Complete accessories included*

Technical Data:

<i>Lamps</i>	<i>2 x Halogen 250 W</i>
<i>Color temperature</i>	<i>approx. 3200 K</i>
<i>Light intensity control</i>	<i>mechanical aperture (infinitely variable)</i>
<i>Line voltage</i>	<i>115/230/240 V AC</i>
<i>Frequency</i>	<i>50/60 Hz</i>
<i>Power consumption</i>	<i>260 W</i>
<i>Operating temperature</i>	<i>+10°C to +40°C</i>
<i>Dimensions (WxHxD)</i>	<i>300 x 145 x 320 mm</i>
<i>Weight</i>	<i>approx. 10 kg</i>

Subject to technical changes.