

62


Galle, Niere
Gall Bladder, Kidney
Vésicule biliaire, Rein
Vesicula biliar, Riñon
Vesicola biliare, Rene




62.101.23
Fig.1




62.103.23
Fig.1



cm DESJARDINS
 24,0 **62.103.24**
 22,5 **62.103.25**



62.103.24



62.103.25

sehr schlankes Modell
 very slender pattern
 modèle très fin
 modelo muy delgado
 modello molto delicato

cm DESJARDINS
 23,0 **62.101.23**
 23,0 **62.103.23**


Gallensteinzangen
Gall Stone Forceps
Pincés pour calculs biliaires
Pinzas para cálculos biliares
Pinze per calcoli biliari

dimeda[®]

SURGICAL INSTRUMENTS



cm CZERNY
25,0 62.105.25



cm BLAKE
21,0 62.107.20
23,0 62.107.23

cm BLAKE
21,0 62.106.20

Gallensteinzangen
Gall Stone Forceps
Pincés pour calculs biliaires
Pinzas para cálculos biliares
Pinze per calcoli biliari

dimeda[®]

SURGICAL INSTRUMENTS



cm MAYO BLAKE
20,0 62.110.20



cm MAYO BLAKE
20,0 62.111.20



cm MIXTER
22,0 62.117.22

Gallenkanaldilatatoren
 Gall Duct Dilators
 Dilatateurs pour les voies biliaires
 Dilatadores para conducto biliar
 Dilatatori per condotti biliari

dimeda[®]


SURGICAL INSTRUMENTS

Gallensteinlöffel
 Gall Stone Scoops
 Curettes pour calculs biliaires
 Cucharillas para coger cálculos biliares
 Cucchiai per calcoli biliari



cm MOORE
 28,0 **62.125.28**

federnd
 elastic
 élastiques
 elásticas
 elastiche



62.130.13


62.130.18

62.130.21

cm DESJARDINS
 28,0 **62.130.13**
 28,0 **62.130.18**
 28,0 **62.130.21**

Gallensteinsonde
 Gall Stone Probes
 Sondes à calculs biliaires
 Sondas-exploradoras para cálculos biliares
 Sonde per calcoli biliari

biegsam
 malleable
 souples
 maleables
 malleabili



cm MOYNIHAN
 34,0 **62.142.34**

biegsamer Schaft
 malleable shaft
 tige malléable
 vástago maleable
 gambo flessibile



BAKES



62.150.01
 32 cm
 Fig. 1
 Ø 1 mm



62.150.02
 32 cm
 Fig. 2
 Ø 2 mm



62.150.03
 32 cm
 Fig. 3
 Ø 3 mm



62.150.04
 32 cm
 Fig. 4
 Ø 4 mm



62.150.05
 32 cm
 Fig. 5
 Ø 5 mm



62.150.06
 32 cm
 Fig. 6
 Ø 6 mm



62.150.07
 32 cm
 Fig. 7
 Ø 7 mm



62.150.08
 32 cm
 Fig. 8
 Ø 8 mm



62.150.09
 32 cm
 Fig. 9
 Ø 9 mm




62.150.10
 32 cm
 Fig. 10
 Ø 10 mm



62.150.11
 32 cm
 Fig. 11
 Ø 11 mm



62.150.12
 32 cm
 Fig. 12
 Ø 12 mm



62.150.13
 32 cm
 Fig. 13
 Ø 13 mm

BAKES

62.150.00


Satz à 9 Stück in Leinentuchtaschen
 set of 9 pieces in canvas pocket
 jeu de 9 pièces dans une poche en toile
 juego de 9 unidades en bolsa de tela
 set di 9 pezzi

62.150.03 - 62.150.11

Gallensteinlöffel
Gall Stone Scoops
Curettes pour calculs biliaries
Cucharas para cálculos biliaries
Pinze per calcoli biliari

dimeda[®]

SURGICAL INSTRUMENTS



cm FERGUSSON
21,5 62.170.21



cm MAYO
28,0 62.175.27



cm MAYO
27,0 62.178.01
27,0 62.178.02
27,0 62.178.03

Gallensteinlöffel
Gall Stone Scoops
Curettes pour calculs biliares
Cucharas para cálculos biliares
Pinze per calcoli biliari

dimeda[®]

SURGICAL INSTRUMENTS



federnd
elastic
élastiques
elásticas
elastiche



62.180.17



62.180.21



62.180.24



62.180.27

cm	DESJARDINS
30,0	62.180.17
30,0	62.180.21
30,0	62.180.24
30,0	62.180.27

biegsam
 malleable
 souples
 maleables
 malleabili



62.190.01
 32 cm
 Fig. 000



62.190.02
 32 cm
 Fig. 00



62.190.03
 32 cm
 Fig. 0



62.190.04
 32 cm
 Fig. 1



62.190.05
 32 cm
 Fig. 2



62.190.06
 32 cm
 Fig. 3



62.190.07
 32 cm
 Fig. 4



62.190.08
 32 cm
 Fig. 5

Fig	mm	LUER-KOERTE
000	Ø 2,5	62.190.01
00	Ø 3,2	62.190.02
0	Ø 4,2	62.190.03
1	Ø 5,5	62.190.04
2	Ø 6,7	62.190.05
3	Ø 8,0	62.190.06
4	Ø 11,0	62.190.07
5	Ø 15,0	62.190.08

Ureterklemmen
Ureter Clamps
Pinces-clapms uretère
Pinzas clamps para ureter
Pinze per uretere

dimeda[®]


SURGICAL INSTRUMENTS




cm MAYO
32,0 62.210.32




cm LURZ
21,0 62.221.21




cm LURZ
21,0 62.223.21




62.231.01
 23 cm
 Fig. 1




62.231.02
 21 cm
 Fig. 2




62.231.03
 19 cm
 Fig. 3



62.231.04
 18,5 cm
 Fig. 4




62.231.05
 19 cm
 Fig. 5




62.231.06
 19 cm
 Fig. 6

Fig.	RANDALL
1	62.231.01
2	62.231.02
3	62.231.03
4	62.231.04
5	62.231.05
6	62.231.06

feines Modell
delicate pattern
modèle délicat
modelo fino
modello delicato



62.231.23
23 cm
Fig. 1




62.231.24
24,5 cm
Fig. 2

Fig.	RANDALL
1	62.231.23
2	62.231.24

Nierensteinzangen
 Kidney Stone Forceps
 Pincés pour les calculs du bassinét
 Pinzas para cálculos renales
 Pinze per calcoli renali


dimeda[®]

SURGICAL INSTRUMENTS




62.233.16
 16 cm
 Fig. 1

Nierenfistelzangen
 Nephrostomy Forceps
 Pincés à néphrostomie
 Pinzas para nephrostomie
 Pinze per nefrostomia



62.235.01
 21 cm
 Fig. 1

für Kinder
 for children
 pour enfants
 para niños
 per bambini



62.233.17
 17 cm
 Fig. 2



62.235.02
 21 cm
 Fig. 2

Fig. HEIMING
 1 **62.233.16**
 2 **62.233.17**

Fig. ELSAESSER
 1 **62.235.01**
 2 **62.235.02**

dimeda®

SURGICAL INSTRUMENTS